

Follow these instructions to properly position your Serenity headgear and mask. The device must be turned on in order for the mask to seal to face properly. There should be minimal or no leaks. **NOTE:** *Refer to instruction manual for more detailed instructions.*

STEP 1




FIG A: Pull the assembled mask and headgear over your head.

FIG B: Place the mask cushion over your nose, and ensure that the bottom straps are below your ear. Attach Quick Release buckle.

STEP 2





FIG A: Check the position of the mask on your face and adjust headgear tension starting with the Bottom Left (BL) and Bottom Right (BR) straps.

FIG B: Adjust the Top Left (TL) and Top Right (TR) straps until the forehead pads rest comfortably against your face.

STEP 3



To adjust the position of the forehead rest, simultaneously press the tabs inward on the forehead rest and pull the forehead rest upward away for the mask body.

STEP 4



PUSH the forehead rest down in the best location and lock into place. Put the mask in position over your nose and ensure it is parallel to your face.

STEP 5



Attach CPAP tubing to 15mm tubing w/swivel adapter on Serenity mask.

STEP 6



Place mask on your face and turn on your CPAP device at the prescribed pressure to check for mask leaks.