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
1 Introduction

1.1 Intended use

The nasal masks are used for treating sleep apnea and for non-invasive and non-life-sustaining ventilation of patients with respiratory insufficiency. They serve as a connecting element between the patient and the therapy device.

1.2 Contraindications

In the case that any of the following symptoms are present the mask may not be used, or only with particular caution:

erosions and ulcerations, skin allergies, flush of the facial skin, pressure marks on the face, claustrophobia, anxiety, deformities of the face or nasopharynx, intake of medication which may induce vomiting, the necessity for immediate intubation. Please also observe the contraindications in the instructions for use of your therapy device.

1.3 Side effects

Nasal congestion, dry nose, morning xerostomia (dry mouth), feeling of pressure in the sinuses, irritated conjunctiva, skin rashes, pressure marks on the face, irritating noises when breathing.

2 Safety

Please read these instructions carefully. They are a constituent part of the mask and must be available at all times.

Risk of injury from the supply of oxygen!

Oxygen can deposit in clothing, bed linen, and hair. In conjunction with smoking, naked flames, or electrical equipment, it can cause fires and explosions.

- ⇒ Do not smoke.
- ⇒ Keep away from naked flames.

Risk of injury from re-inhaling CO₂!

If the mask is used incorrectly, CO₂ may be re-inhaled.

- ⇒ Only apply the mask for a prolonged period if the therapy device is running.
- ⇒ Only use the mask within the quoted therapy pressure range.
- ⇒ Do not use the mask on patients who cannot take it off themselves.

Risk of injury from slipping of the mask!

- ⇒ If the mask slips or falls off, the therapy is ineffective.
- ⇒ Monitor patients with restricted spontaneous respiration.
- ⇒ Activate low pressure/leakage alarms on the therapy device.

Risk of injury from anesthetic gases!

Anesthetic gases can escape through the exhalation system and put third parties at risk.

- ⇒ Do not use the mask during anesthesia.

3 Product description

3.1 Overview

Illustrations of the individual components can be found on the front cover.

Key

1. Headgear
2. Forehead cushion
3. Forehead support
4. Retaining ring
5. Elbow
6. Plug for pressure measurement port
7. Pressure measurement port
8. Rotating sleeve
9. Mask body
10. Headgear clip
11. Mask cushion

3.2 Compatible devices

With some combinations of devices, the actual pressure in the mask will not correspond to the prescribed therapy pressure, even if the therapy device is displaying the correct pressure. Have the combination of devices adjusted by a physician or specialist dealer so that the actual pressure in the mask corresponds to therapy pressure.

3.3 Exhalation system

The mask has an integrated exhalation system. Retaining ring and mask body are shaped in such a way that a gap forms between these parts. The exhaled air can escape through this gap.

3.4 Pressure measurement port

The mask has a pressure measurement port which can be used for measuring therapy pressure or supplying oxygen. If you are not using this, close it with the plug for pressure measurement port.

4 Operation

The following figures show you how to adjust, apply, remove, dismantle and assemble the mask:

Figure 1: Putting on the mask

Figure 2: Adjusting the mask

Figure 3: Removing the mask

Figure 4: Dismantling the mask

Figure 5: Assembling the mask

5 Hygienic preparation

⚠ WARNING

Risk of injury from inadequate cleaning!

Residues may congest the mask, impair the integrated exhalation system, and jeopardize therapy success.

⇒ For patients with a compromised immune system or particular background of illness, disinfect mask components daily following consultation with the physician.

5.1 Cleaning the mask

1. Dismantle the mask (see figure 4).
2. Clean mask in accordance with the table below:

Action	Daily	Weekly
Wash mask components with warm water and a mild detergent.	X	
When washing mask parts (exception: mask cushion and forehead cushion), clean thoroughly with a cloth or a soft brush.		X
Wash headgear by hand.		X

3. Rinse all parts with clear water.

4. Allow all parts to air-dry.

5. Perform a visual inspection.

6. If necessary: Replace damaged parts.

7. Assemble the mask (see figure 5).

i

Discolorations of mask components do not impair the functionality of the mask.

5.2 Change of patient

You must subject the mask to a hygienic preparation in the event of a change of patient. For information on hygienic preparation in the event of a change of patient, please refer to our Internet brochure at www.weinmann-medical.com. If you do not have Internet access, please contact Weinmann.

6 Disposal

You can dispose of all parts in domestic waste.

7 Faults

Fault	Cause	Remedy
Facial pressure pain.	Mask fits too tightly.	Loosen headgear.
Draft in the eyes.	Mask fits too loosely.	Tighten headgear.
	Mask does not fit.	Contact your authorized dealer.
	Mask is not adjusted correctly.	Re-adjust mask (see figure 2).
	Mask cushion is damaged.	Replace mask cushion.
Therapy pressure is not reached.	Tube system is leaking.	Check connection and fit of tubes.
	Air is escaping from pressure measurement port.	Close pressure measurement port with plug (included in scope of supply).

8 Technical data

Product class according to 93/42/EEC	IIa	
Dimensions (W x H x D)	80 mm x 130 mm x 90 mm	
Weight	72 g	
Therapy pressure range	4 hPa - 30 hPa	
Hose connection:	Ø 22 mm (male)	
Cone according to EN ISO 5356-1		
Pressure measurement port	Ø 4 mm	
Width of headgear clips	20 mm	
Temperature range:		
Operation		+5°C to + 40°C
Storage		-20°C to + 70°C
Flow resistance at 50 l/min		0.16 hPa
Flow resistance at 100 l/min		0.64 hPa

Quoted two-figure noise emission value according to ISO 4871:	15 dB(A)
- Sound pressure level	23 dB(A)
- Uncertainty factor	3 dB(A)
Service life	Up to 12 months ¹
Standards applied	EN ISO 17510-2:2009


¹The materials of the mask deteriorate if exposed to e.g., aggressive detergents. In individual cases, therefore, it may be necessary to replace mask parts sooner.

CE 0197

The right to make design modifications is reserved.

9 Pressure/flow characteristic curve

The pressure/flow characteristic curve shows the leak flow as a function of therapy pressure.

10 Materials

Mask component	Material
Mask cushion, forehead cushion	Silicone, silicone coating
Plug for pressure measurement port, spring of forehead support	Silicone
Retaining ring (hard component)	PP (polypropylene)
Retaining ring (soft component)	TPE (thermoplastic elastomer)

12 Declaration of conformity

Weinmann Geräte für Medizin GmbH + Co. KG (Kronsauweg 40, 22525 Hamburg, Germany) hereby declares that the product complies with the relevant regulations of Directive 93/42/EEC governing medical devices. The complete text of the declaration of conformity can be found at: www.weinmann-medical.com

Mask component	Material
Headgear clip	PA (polyamide), POM (polyoxymethylene)
Forehead support, rotating sleeve, mask body, elbow	PA (polyamide), PC (polycarbonate)
Headgear	Elastane, polyester, PU (polyurethane), UBL loop, cotton, PA (polyamide)

All mask components are free from latex, PVC (polyvinyl chloride) and DEHP (diethylhexyl phthalate).

11 Warranty

Weinmann gives the customer a limited manufacturer warranty on new genuine Weinmann products and any replacement part fitted by Weinmann in accordance with the warranty conditions applicable to the product in question and in accordance with the warranty periods from date of purchase as listed below. The warranty conditions can be downloaded from www.weinmann-medical.com on the Internet. We can also send you the warranty conditions on request. Please note that any claim under warranty and liability will be void if neither the accessories nor original replacement parts recommended in the instructions for use are used.

In the event of a claim under warranty, please contact your specialist dealer.

Product	Period of guarantee
Masks including accessories	6 months

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