

VECTOR *et* BILEVEL



VECTOR *et* BILEVEL

The device is a therapy device in the BILEVEL version for treating obstructive sleep apnea.

Trigger device

The respiration trigger detects the patients efforts to breath in and out and signals these efforts to the control device. A flow-based trigger is integrated which can be set separately for inspiration and expiration. While the pressure increases during inspiration it decreases during expiration. In connection with an optimum flank adjustment this results in a ventilation which is very comfortable for the patient and which gives the patient the impression as if the device is following the natural respiration without delay.

To generate the necessary pressure, the device is provided with an electronically controlled blower. In order to keep the stress for the patient as low as possible, the blower has been designed with a high reserve capacity and a quick control response.

In addition to the mask alarm, the device offers another safety function, the power failure alarm. It automatically restarts the motor after power restoration.

For added comfort an alarm clock is integrated in the device. Further comfort functions include a soft start ramp, the automatic start-stop function and the adjustable display brightness.

Technical Data

Dimensions (w x d x h)	230 x 212 x 107 mm
Mass	approx. 2,0 kg
Power supply	100 V ~...240 V ~, 50...60 Hz 24 V DC / 2,1 A
Pressure range	4 ... 20 hPa (mbar)
Max. limit pressure in the event of an error	< 40 hPa (konstruktiv bedingt)
Power consumption (24 V DC)	Maximum 50 W Standby < 5 W at 8 hPa (without humidifier) 10 W at 28 hPa (with humidifier, level 5) 30 W
Power consumption (power supply)	Maximum 50 W Standby < 5 W at 10 hPa (without humidifier) 12 W at 32 hPa (with humidifier, level 5) 33 W
Short time pressure variance at 4 hPa	10 bpm = ± 0,10 hPa 15 bpm = ± 0,19 hPa 20 bpm = ± 0,28 hPa
Short time pressure variance at 8 hPa	10 bpm = ± 0,10 hPa

	15 bpm = ± 0,19 hPa
	20 bpm = ± 0,28 hPa
Short time pressure variance at 12 hPa	10 bpm = ± 0,10 hPa
	15 bpm = ± 0,19 hPa
	20 bpm = ± 0,28 hPa
Short time pressure variance at 16 hPa	10 bpm = ± 0,10 hPa
	15 bpm = ± 0,19 hPa
	20 bpm = ± 0,28 hPa
Short time pressure variance at 20 hPa	10 bpm = ± 0,10 hPa
	15 bpm = ± 0,19 hPa
	20 bpm = ± 0,28 hPa
Values determined with and without AquaTRENDuni	
Long time pressure variance	< 0,3 hPa
Display accuracy pressure	0,5 hPa
Average noise level	< 25 dB(A) bei 10 hPa Betrieb in 1 m Abstand (entspricht einem Schalleistungspegel von 33 dB(A))
Air delivery volume	> 170 l/min at 4 hPa therapy pressure > 160 l/min at 8 hPa therapy pressure > 150 l/min at 12 hPa therapy pressure > 130 l/min at 16 hPa therapy pressure > 120 l/min at 20 hPa therapy pressure
Operating temperature	+ 5 °C ... + 40 °C
Storage temperature	- 20 °C ... + 70 °C
Air humidity	15 ... 95 % rel. air humidity (operation and storage)
Operating altitude	- 400 m ... 3500 m (1060 hPa ... 700 hPa)
Filtration efficiency up to 1 µm	99,5 %
CE-mark	according to EC Directive 93/42/EEC
Product class according to the MDD	Ila

The technical specifications are subject to technical modification.