N4 Nasal Mask User Manual

N4 Nasal Mask

Thank you for choosing BMC's N4 Nasal Mask. The Nasal Mask is comfortable, flexible and quiet. It is designed to minimize contact with your face, thus ensuring that you feel comfortable during therapy. This user manual provides you with the information you need for the correct use of your mask.

Intended Use

The N4 Nasal Mask channels airflow non-invasively to a patient from a positive airway pressure device such as a continuous positive airway pressure (CPAP) or Bi-level system.

The N4 Nasal Mask is:

- To be used by adult patients (> 66 lb / 30 kg) for whom positive airway pressure therapy has been prescribed.
- Intended for single-patient re-use in the home environment and multi-patient re-use in the hospital / institutional environment.


CAUTION: In the US, Federal law restricts the sale of this device to sale by or on the order of a physician.

Before Using the Mask

⚠️ WARNINGS

- The vent holes must be kept clear.
- This mask should only be used with CPAP or bilevel devices recommended by a physician or respiratory therapist.
- The mask should not be used unless the device is turned on. Once the mask is fitted, ensure the device is blowing air.
- Explanation:** CPAP and bilevel devices are intended to be used with special masks (or connectors) which have vent holes to allow continuous flow of air out of the mask. When the device is turned on and functioning properly, new air from the device flushes the exhaled air out through the mask vent holes. However, when the device is not operating, insufficient fresh air will be provided through the mask, and the exhaled air may be rebreathed. Rebreathing of exhaled air for longer than several minutes can, in some circumstances, lead to suffocation. This applies to most models of CPAP or bilevel devices.
- This mask should not be used on patients who are uncooperative, obtunded, unresponsive, or unable to remove the mask.
- Follow all precautions when using supplemental oxygen.
- Not applicable for silicone-allergy user.
- Oxygen flow must be turned off when the CPAP or bilevel device is not operating, so that unused oxygen does not accumulate within the device enclosure and create a risk of fire.
- Oxygen supports combustion. Oxygen must not be used while smoking or in the presence of an open flame. Only use oxygen in well ventilated rooms.
- At a fixed flow rate of supplemental oxygen flow, the inhaled oxygen concentration varies, depending on the pressure settings, patient breathing pattern, mask, point of application and leak rate.
- The technical specifications of the mask are provided for your clinician to check if it is compatible with the flow generator. If it is used beyond technical specifications or used with incompatible devices, the seal and comfort of the mask may not be effective, optimum therapy may not be achieved, and leak, or variation in the rate of leak, may affect the function of the flow generator.
- Stop using the N4 Nasal Mask and consult your physician or sleep therapist, if you have ANY adverse reaction to the use of the mask.
- Refer to your CPAP or bilevel device manual for details on settings and operational information.
- Remove all packaging before using the mask.

Components of the Mask


Fitting the Mask

- Attach the headgear straps by threading them through the gaps of forehead frame and clips.
- Hold the mask over your nose and pull the headgear over your head.
- Bringing the lower headgear straps below your ears, connect the clips to the Frame assembly.
- Adjust the upper headgear straps evenly until they are comfortable and firm, and reattach the mask.
- Adjust the lower headgear straps evenly until they are comfortable and firm, and reattach the mask.

Disassembling the Mask

Note:

The Frame assembly cannot be disintegrated.

Cleaning the Mask at Home

Notes:

- The mask and headgear can only be cleaned through hand washing.
- The headgear may be washed without removing the clips.
- When washing the mask, the Elbow, Swivel, Forehead frame and Cushion frame cannot be disassembled from the Frame assembly.

⚠️ CAUTION

- Always follow cleaning instructions and do not use solutions containing bleach, chlorine, alcohol, aromatics, moisturizers, antibacterial agents, or scented oils to clean any part of the system or air tubing, otherwise may damage the mask.
- Exposing any part of the system or tubing to direct sunlight or heat may cause deterioration.
- If any visible deterioration of a component is apparent (cracking, crazing, tears etc), the component should be discarded and replaced.

Daily / After Each Use

- Disassemble the mask components according to the disassembly instructions.
- Thoroughly clean the separated mask components (excluding headgear), by gently rubbing in warm (approx. 30°C) water using mild soap.
- Use a soft bristle brush to clean the vent.
- Rinse all components well with drinking quality water and allow to air dry out of direct sunlight.
- When all components are dry, reassemble according to the reassembly instructions.


Weekly

Handwash the headgear and all components in warm (approx. 30°C) and soapy water. Use pure soap. Rinse the components well and allow them to air dry out of direct sunlight.

Reprocessing the Mask between Patients

This mask should be reprocessed when used between patients. Please refer to the accompanying "Cleaning and Disinfection Instructions" or BMC website www.bmc-medical.com for "Cleaning and Disinfection Instructions". If you do not have internet access, please contact your BMC representative.


Reassembling the Mask

Troubleshooting

Problem	Possible Reason	Possible Solution
Mask won't seal properly or is uncomfortable.	Mask may have been fitted incorrectly or is incorrectly adjusted.	Carefully follow instructions in "Fitting the Mask" section. Make sure the headgear is not over-tightened.
Mask size is wrong.		Consult your clinician.
Mask is not tightened correctly.		Tighten both upper and lower headgear straps.
Cushion is misaligned or creased on face.		Ensure cushion is properly aligned according to the instructions in "Reassembling the Mask". Make sure cushion has no creases in it.
Mask cushion may be dirty.		Clean cushion according to instructions in "Cleaning the Mask at Home".
Mask is too noisy.	Mask is assembled incorrectly.	Reassemble mask.
	Vents are blocked or partially blocked.	Clean the vents according to instructions in "Cleaning the Mask at Home" section.

Technical Specifications

Dead Space Information	Dead space is the empty volume of the mask up to the swivel. The dead space of the mask varies according to cushion sizes but is less than 92 mL.
Therapy Pressure	4 to 25 cmH ₂ O.
Resistance	Drop in Pressure measured (average for 3 sizes) at 50 L/min: 0.2 cmH ₂ O at 100 L/min: 0.8 cmH ₂ O
Sound	DECLARED DUAL-NUMBER NOISE EMISSION VALUES in accordance with ISO 4871. The A-weighted sound power level of the mask is 28 dBA, with uncertainty 3 dBA. The A-weighted sound pressure level of the mask at a distance of 1 m is 20 dBA, with uncertainty 3 dBA.
Environmental Conditions	Operating temperature: +5°C to +40°C (41°F to 104°F) Operating humidity: 0% ~ 93% relative humidity non-condensing Storage and transport: -20°C to +55°C (-4°F to 131°F) Storage and transport humidity: 0% ~ 93% relative humidity non-condensing

Storage

Ensure that the mask is thoroughly clean and dry before storing it for any length of time. Store the mask in a dry place out of direct sunlight.

Disposal

The N4 Nasal Mask does not contain any hazardous substances and may be disposed of with your normal household refuse.

Symbols

System and Packaging

	Caution, consult accompanying documents. Indicates a Warning or Caution and alerts you to a possible injury or explains special measures for the safe and effective use of the device.
	Lot number
	Catalog number
	Not made with natural rubber latex
	Temperature limitation
	Humidity limitation
	Manufacturer
	Authorized Representative in the European Community

Limited Warranty

It is warranted that N4 Nasal Mask, including its Frame assembly, Cushion, Clips and Headgear shall be free from defects in material and workmanship for a period of ninety (90) days from the date of purchase by the initial consumer. To exercise the rights under this warranty, contact the local authorized dealers.

The expected service life of N4 Nasal Mask is one year.

MANUFACTURER:

BMC MEDICAL CO., LTD.
Room 110 Tower A Fengyu Building, No. 115 Fucheng Road, Haidian, 100036 Beijing, PEOPLE'S REPUBLIC OF CHINA

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AB YETKİLİ TEMSİLCİSİ:
Shanghai International Holding Corp. GmbH (Europe)