

Nasal-Aire II Petite

Nasal-Aire® II Petite Interface

A Revolutionary Interface for Pediatric and Smaller Adult Patients

As one of the only interfaces on the market **approved for use by pediatric patients**, the Nasal-Aire® II Petite with lightweight Micro™ Headgear fits comfortably and easily on even the smallest faces.

Small and quiet, the Nasal-Aire II Petite nasal cannula mask was designed for patients over seven years of age who require positive airway pressure therapy for Sleep Disordered Breathing.

Sharing the same comfort and compliance benefits of the Nasal-Aire technology, its simplistic, easy-to-use design makes it ideal for younger patients.

Performance Features

- Virtually leak-free and ergonomically designed to fit the smallest of faces
- Minimalistic design virtually eliminates
 pressure points on the face allowing the
 interface to be worn for extended periods
 of time without discomfort
- Nasal cannula design optimizes seal within the nose, offers an unobstructed line of sight and helps to reduce feelings of claustrophobia
- Headgear with secure fit design allows the user to sleep in most any position comfortably
- **Built-in swivel** prevents tubing from pulling and kinking, allowing a more comfortable, uninterrupted night's sleep

Design Features

- Designed for patients over 7 years of age or weighing more than 40 lbs.
- Hanging tubing reduces moisture build-up underneath the nose
- Lack of nasal bridge and forehead contact virtually eliminates sores and air leaks near eyes
- Easy to fit, ergonomically designed Micro headgear with patented elbow feature reduces facial pressure points for added comfort

- Patients can rest comfortably in various positions without the constraint of traditional nasal masks and bulky headgear
- Simplistic design requires less parts than traditional nasal mask systems and makes interface easy for pediatric patients to assemble and clean
- Available in five sizes smaller than the Nasal-Aire II extra small
- 100% latex free soft, pliable silicone

Reimbursement Information

PRODUCT DESCRIPTION	CODE DESCRIPTION	HCPCS CODING	FREQUENCY	ALLOWABLE RANGE
Nasal-Aire II Petite	Nasal interface (mask or cannula type) Used with PAP device,with or without head strap	A7034	1 per 3 months	\$92.57 - \$108.90
Nasal-Aire II Petite Insert	Replacement pillows for nasal application device, pair	A7033	2 per 1 month	\$22.36 - \$26.30
Headgear	Headgear used with PAP device, each	A7035	1 per 6 months	\$31.27 - \$36.79

The noted allowable rates are based off of the January 2012 HCPCS Fee Schedule. All information presented is for reference purposes only.

Accessories and Replacement Parts

PRODUCT DESCRIPTION	PART NUMBER
Nasal-Aire II Petite, Size A Nasal-Aire II Petite, Size B Nasal-Aire II Petite, Size C Nasal-Aire II Petite, Size D Nasal-Aire II Petite, Size E Nasal-Aire II Petite – Headgear Nasal-Aire II Petite – Cannula Only, Size A Nasal-Aire II Petite – Cannula Only, Size B Nasal-Aire II Petite – Cannula Only, Size C Nasal-Aire II Petite – Cannula Only, Size D Nasal-Aire II Petite – Cannula Only, Size E	PA401 PB402 PC403 PD404 PE405A P430 PA411 PB412 PC413 PD414 PE415A
Sizing Gauge – Nasal-Aire Petite	X4007

Nasal-Aire II Petite Cannula Measurements

SIZE	MAJOR OD (MM)	MINOR OD (MM)
А	13.61	10.37
В	12.57	9.57
С	11.48	8.97
D	10.78	8.33
Е	10.03	7.88

Use of a sizing gauge is recommended to achieve the best fit for the patient.



DeVilbiss Healthcare LLC | 100 DeVilbiss Drive | Somerset, PA 15501 USA 800.338.1988 | 814.443.4881 | www.DeVilbissHealthcare.com