



Respiratory Therapy

XT Heated Humidifier

Deliver a Proper Care

XT series provides an add-on option for patients who require heated humidity with their CPAP therapy for increasing comfort of therapy.

Small and compact design that work with our new XT series family. It allows patients to choose the most comfortable level of humidification with the rotary knob. The power supply is powered by the unit with jumper cord, so no separate power cord is required. The platform that is perfect for stacking with CPAP to stay in place. All in all, the compact, lightweight, streamlined design of XT heated humidifier work well wherever you want it to be.





Features:

- Integrated with XT series CPAP system.
- Easy to use with the rotary knob for heating level adjustment. .
- The see through reservoir allows for easy monitoring of the water level.
- Water chamber with removable base plate allows easy cleaning
- Safe & reliable with overheating protection mechanism.
- Conspicuous indication to show Power ON and Heater ON.



Specifications:

Model	XT Heated Humidifier
Reservoir Capacity:	450ml
Heater Settings:	1 to 6 from 30 °C to 60°C (86°F to 140°F)
Humidity Output:	10-40 mgH2O/L
Pressure Drop:	0.2 cmH2O @ 60LPM
Dimensions (L x W x H):	10" x 6.8" x 5.8" (25.5 cm x 17.2 cm x 14.6 cm)
Weight:	1.37 lb (620g)
Power Supply:	120VAC, 50/60Hz (For 120V system) 220V-240VAC, 50/60Hz (For 230V system)

Specifications are subject to change for improvement without notice.

ACCESSORIES



XT series CPAP System includes:

- CPAP
- Power cord
- 6ft Hose
- User's manual



6 Ft. CPAP Hose

- 22mm OD endcaps will fit standard CPAP/BiPAP setups from most manufacturers.
- Create a durable and flexible, whisper quiet connection between mask and CPAP machine.



Carrying Bag

- With soft padding provides perfect protection to CPAP device and accessories for traveling
- Inner pocket for placing masks, tubing and power cord arrangement



Easy Compliance 2.1 Data Management Software

- Exclusively for XT-III CPAP
- Daily pressure chart and statistics report of compliance data during a selected period makes patient management easy and efficient
- USB cable include

APEX MEDICAL CORP.

9, Min Sheng St. Tu-Cheng,
 Taipei County, 236, Taiwan, R.O.C.
 TEL: 886-2-22683100
 FAX: 886-2-22686525
 E-mail:marketing@apexmedical.com.tw
 www.apexbrand.com

APEX MEDICAL S.L.

Máximo Aguirre 18 Bis, 8ª planta
 48011 Bilbao (Vizcaya)
 Tel: (+34) 94.470.64.08
 Fax: (+34) 94.470.64.09
 E-mail: info@apexmedical.es
 www.apexbrand.com