Nasal-Aire® II

The Nasal-Aire® II is a nasal cannula style interface for CPAP therapy. This lightweight interface wears like a familiar oxygen cannula. The Nasal-Aire® II All Sizes Kit contains all of the components required to properly size the patient, eliminating the hassle of opening multiple products to get the right fit.

Product Features

- Cannula Style interface with patented sealing within the nose
- Virtually eliminates pressure points on the face and bottom of nose
- Features patented L-Flow[®] technology to deliver flow to patient with less turbulence
- Headgear tensioning not required as cannula creates the seal
- Available in an All Sizes Kit or in individual Size Systems, including a Petite Sized product line





Cannulas available in 6 different sizes. Each size is designed to be 9% larger than the previous size.

(Note that the MP size fills an 18% sizing gap between the MD and LG Cannulas.)

Order Information

| Part No. | Description | HCPCS Code | Code Description | HCPCS Frequency |
|--|--|---------------------|--|------------------------------------|
| K2A | Nasal-Aire II – All Sizes Kit | A7034 + A7035 | Nasal interface used with PAP device; Headgear used with PAP device | 1 per 3 months + (see below) |
| XL306 LG305 MP304 MD303 SM302 XS301 | Nasal-Aire II Individual Systems – XL LG MP (Medium Plus) MD SM XS | A7034 + A7035 | Nasal interface used with PAP device; Headgear used with PAP device | 1 per 3 months + (see below) |
| NAT326 NAT325 NAT324 NAT323 NAT322 NAT321 | Nasal-Aire II Frame Systems – XL LG MP (Medium Plus) MD SM XS | A7034 | Nasal interface used with PAP device (no headgear included) | 1 per 3 months |
| XL06 LG05 MP04 MD03 SM02 XS01 | Replacement Cannulas – XL LG MP (Medium Plus) MD SM XS | A7033 | Cannula for nasal interface | 2 pair per 1 month |
| X2105 | Replacement Headgear | A7035 | Headgear used with PAP device | 1 per 6 months |
| X2007 | Sizing Gauge for Nasal-Aire II | - | - | - |

