

Swing into Quattro FX

Achieving an effective seal is the key to successful therapy; but it needn't be complicated. Quattro FX introduces the two-step Swing Seal fitting technique: for a secure seal first time, every time.

- Firstly position the cushion at an angle on the nasal bridge ...
- ... then **swing** the lower half of the mask down against the face to **seal**.



Simply swing, and seal.

RESMED

Less mask. More you.

www.quattrofx.com





To get the best out of your therapy, pair your Quattro FX with a sleep therapy system from ResMed's S9TM Series.

Visit www.resmed.com for more information.

Manufacturer:ResMed Pty Ltd 1 Elizabeth Macarthur Drive Bella Vista NSW 2153 Australia. Distributed by: ResMed Corp 9001 Spectrum Center Boulevard San Diego CA 92123 USA. ResMed (UK) Ltd 96 Milton Park Abingdon Oxfordshire 0X14 4RY UK. See www.resmed.com for other ResMed locations worldwide. S9, Quattro and Spring Air are trademarks of ResMed Pty Ltd and are registered in U.S. Patent and Trademark Office. © 2010 ResMed. Specifications may change without notice. 1014221/1 2010-08

RESMED

Less mask. More you.

User Information





The unobtrusive alternative

It may be called a full face mask, but it certainly doesn't look like one! Compared to conventional full face masks, Quattro FX is light and unobtrusive, giving you an effective alternative, and a sense of freedom, in full face therapy.

Support and freedom

Quattro FX provides the perfect balance of stability and flexibility. The contoured headgear design stabilises the mask on your face, and the Spring Air™ cushion is soft, comfortable and secures your seal even if you move around in your sleep.

A mask you'll feel good about

Being on therapy doesn't mean losing your sense of self. Even though it's a full face mask, Quattro FX starts mid-way on the nasal bridge, giving you a clear line of sight and a feeling of openness. So you'll enjoy a good night's sleep and wake up to a new day with a renewed outlook on your mask and your health.

Quattro[™] FX
FULL FACE MASK