

The compact full face mask designed to easily fit everyone

The AirFit F20 features an innovative InfinitySeal[™] silicone cushion: tested to perform at high therapy pressures, it's designed to comfortably fit patients of all facial shapes and sizes. It will appeal to patients looking for a less claustrophobic mask that delivers a clear line of sight.



In ResMed international studies, AirFit F20 fit 96.5% of all patients tested.^{1, 2, 3, 4} This aims to make it easier for you to set patients up right the first time and send them home with confidence.

Robust seal

This mask features an adaptive cushion that has been tested to prevent blow-outs, accommodate patient movement and provide a robust seal across high pressures.



Comfortable to use

The flexible fabric-lined frame adds further comfort while delivering all-round stability and greater visual freedom.

With the InfinitySeal silicone cushion, the AirFit F20 is designed to adapt to the unique facial contours of each patient, so you can be confident of providing a secure and comfortable mask fit.



Visual freedom

Let patients see and do more in bed with an under-the-eye frame that combines all-round stability and less facial contact.

Soft and gentle

The plush headgear is specifically designed for extra softness and to help ensure patients settle into therapy quickly and comfortably.

Fast fitting time

Less than one minute is needed to fit AirFit F20.¹ Taking the mask on and off is convenient and intuitive, with magnetic clips that guide the headgear to the frame in seconds.

Dependable seal

Set patients up with a reliable and effective seal thanks to the InfinitySeal silicone cushion – designed to prevent blow-outs, tolerate misfitting and seal well under variable pressures.

Soft and flexible frame

Designed to support the widest range of patients, the soft, flexible frame adapts to different face shapes while its integrated padding adds comfort.

Flexible connection and movement

The elbow's quick release feature enables easy connection or disconnection, while its 360° rotation capability allows for greater freedom of movement.






Designed specifically for women, it offers a smaller size to suit female facial contours.



In ResMed international studies, AirFit F20 was found to fit 96.5% of all patients tested.^{1, 2, 3, 4}



Minimal parts make the AirFit F20 easy to use and clean




Sizing tips

AirFit F20 has three interchangeable cushion sizes - small, medium and large. The modular design of the AirFit F20 means each size cushion easily clicks into the same mask frame, so you can swap between them with ease for more flexibility.

We recommend that you refer to the mask fitting template to gauge which size is the most suitable for each of your patients.

You can download it on ResMed.com/AirFitF20




Product codes

	APAC	ANZ	China	Japan
AirFit F20: Small	63420	63430	63435	63425
AirFit F20: Medium	63421	63431	63436	63426
AirFit F20: Large	63422	63432	63437	63427
AirFit F20 for Her: Small	63423	63433	63438	63428
AirFit F20 for Her: Medium	63424	63434	63439	63429

APAC: English, Simplified Chinese, Traditional Chinese



86% of patients preferred AirFit F20's seal to that of the market-leading full face mask.¹



Evolving cushion technology: variable thickness designed for more comfort and stability

