

family cares



Specification

Size	101*64*44 mm
Weight	145g
Air flow	3.5 CFM
Power consumption	5 W
Power supply	DC 4.5V ~ 6.5V
Storage temperature	-20 °C ~ +55 °C
Operation temperature	-5 °C ~ +45 °C

ResFit Disinfector

Be Focused on Efficient Disinfection.



The disinfector is designed to kill the germs efficiently for ensuring an unpolluted usage.

BMC Medical Co., Ltd.

Room 110, Tower A Fengyu Building, No.115 Fucheng Road, Haidian,
100036 Beijing, PEOPLE'S REPUBLIC OF CHINA
Tel: +86-10-5166 3880 Fax: +86-10-5166 3880 ext. 810 en.bmc-medical.com
SC-20-02-03-31-001-Ver1.0

Technical principle


Ozone is recognized as a powerful oxidant and is used as common disinfectant, which is applicable for the necessities and household environment.
 The disinfectant releases highly concentrated ozone by electrolyzing the sourced natural air.





Product display




Key Features


 99.99%
 Disinfection rate


 Cute body &
 Easy operation


 USB power supply


 2 modes available

Using Guide



Product component

(Note: No battery in the device)

