
Osleepnet[®] Ascend[®] Nasal replaceable AIRgel cushion masks

featuring Advanced Cushion Technology

Rise to a higher level of comfort and convenience with The Ascend nasal mask: our first-ever nasal mask featuring a replaceable cushion. Our patented AIR^ogel replaceable cushion is easily detached and reattached for easy cleaning and cushion replacement, saving both time and money.


Ascend to a Higher Level of Comfort and Convenience



What makes our Ascend Nasal Mask different?

- Our easy-to-use replaceable cushion allows you to extend the life of your mask.
- *Flexible Spacebar* is designed to provide superior stability without contacting the skin.
- Active Headgear Connector moves when the patient moves to help maintain an effective seal.
- Available in three sizes with simple and breathable quick-release EZFit headgear.

Advanced CushionTechnology(A.C.T.)



Easy-to-Replace Cushion





Ascend Nasal Masks with EZFit Disposable Headgear Specifications:

| Product Number | Component | Size | Mask Weight with cushion | Mask Dimensions with cushion |
|-------------------|------------------------------|-----------|--------------------------------|---------------------------------|
| 50174 | Mask System & headgear | S, M, & L | NA | NA |
| 50201 | AIR [•] gel Cushion | Small | 57.9g | 3.70''D × 2.75''W × 5.40''H |
| 50202 | AIR [•] gel Cushion | Medium | 58.9g | 3.80"D × 2.95"W × 5.40"H |
| 50203 | AIR [•] gel Cushion | Large | 63.0g | 3.97''D × 3.30''W × 5.40''H |
| 55051 | EZFit Headgear | NA | NA | NA |





Also available as a Full Face Mask



Contact us today.

To place your order, call an authorized distributor.

Call I.800.742.3646

Visit www.sleepnetmasks.com to learn more.

Sleepnet Corporation, 5 Merrill Industrial Drive, Hampton, New Hampshire 03842 USA



Scan code to view Ascend Nasal Fitting Tips video.