

Accessories Laboratory

TRENDremote

TRENDremote is for the cabled handling of respiration therapy devices from HOFFRICHTER GmbH. With TRENDremote the sleep laboratory staff can change the running therapy without being next to the therapy device. TRENDremote is exclusively used in sleep laboratories and operated by qualified personnel. TRENDremote supports all therapy devices of HOFFRICHTER GmbH.



Dimensions (w x d x h): 115 x 84 x 25 mm
Mass: 134 g
Order number: 00005903

TRENDinterface


TRENDinterface is intended for conversion of digital signals from the respiration therapy devices of HOFFRICHTER GmbH into analogue signals and their transfer to polysomnography systems. Additional the communication between the respiratory therapy device and a PC with evaluation software or a remote control shall be warranted. TRENDinterface is exclusively used in sleep laboratories by medical qualified personnel.



Dimensions (w x d x h): 115 x 84 x 25 mm
Mass: 110 g
Order number: 00005902

DATA box

The DATA box is an accessory device which can be used in connection with the respiratory therapy device series TREND II, VECTOR et and point. It contains a removable memory card and serves to summarise and store the independent transfer of device and therapy data. The memory card can be sent to the sleep laboratory or doctor in charge by the patient, independent of the therapy device. After successful data assessment by medical personnel, the e.g. modified therapy parameters can be copied on the memory card with the help of the DATA box and finally send to the patient again.



Dimensions (w x d x h): 115 x 84 x 25 mm
Mass: 110 g
Order number: 00005982

TRENDset

TRENDset is the communication basis between the Hoffrichter respiration therapy devices and the PC's in sleep laboratories. The data exchange can be done directly via an interface wire or indirect with a SD-Card and the DATA box. This a comfortable possibility for the medical and technical stuff to have a remote control, a therapy analysis and a device check in one software.



Accessories

Accessories Ventilation

Trolley

The trolley is a mobile stand for our ventilator CARAT. The trolley is especial tilt resistant because of the special construction. Optional you can assemble a hinged bracket for the optimal guiding of the therapy tube and a holder for infusion equipment. On the equipment rack is a DIN rail for other accessories. Order number:

Trolley: 00004020
Hinges bracket: 00004021
Tray: 00004022
Infusion pole: 00004023



Alarm box

The remote alarm box is a wire alarm system signalling alarms of the CARAT ventilator for domestic use. Dimensions (w x d x h): 66 x 25 x 125 mm
Mass: 130 g
Order number: 00004035



AKKUPACK CARAT

The external power pack AKKUPACK CARAT is intended as autonomous power supply for the ventilator CARAT. It was especially developed for homecare use, use in hospitals and effort in mobile use. AKKUPACK CARAT enabled an additional using time of 12 hours (in standard adjustments). Dimensions (w x d x h): 220 x 170 x 190 mm
Mass: 6.1 kg
Power supply: PCM120PS24 or SNP-A-129-M
Wide range input: 100 - 250 V ~ at 50 / 60 Hz
Max. power output: 120 W
Order number: 00004030



O₂ - Sensor

The O₂-Sensor measures the oxygen at the air outlet of the ventilator and shows the actual value on the display of the CARAT. Order number: 23000018



Manufacturer:

HOFFRICHTER GmbH
Mettenheimer Straße 12/14
19061 Schwerin
Germany
Phone: +49 385 39925-0
Fax: +49 385 39925-25
E-Mail: info@hoffrichter.de
www.hoffrichter.de

Sales:



Extras make the sound

Accessories

Masks

HOFFRICHTER GmbH offers three different masks in the product portfolio. A respiration therapy nasal mask, a respiration therapy fullface mask and a ventilation mask.

The masks are made of one piece molded silicone rubber. This makes the masks light and comfortable. A standard tubing attachment (22 mm) is turnable in both directions with 360 degrees.

The headgear with Velcro fasteners provide maximum versatility and movement.

The masks are available in three different sizes:

Order number:			
	Small	Medium	Large
Nasal mask	00003440	00003434	00003435
Mouth nasal mask	00003441	00003436	00003437
NIPPV mouth nasal mask	----	00003442	00003438



Nasal mask



Mouth nasal mask



NIPPV mouth nasal mask

Therapy tubes

HOFFRICHTER GmbH offers four different therapy tubes for the different therapy systems. One tube system for the therapy device point without an integrated pressure measuring tube, one tube system for the therapy device series TREND, TREND II, TRENDvent and VECTOR et with integrated pressure measuring tube, an single line patient circuit for CARAT I and II as well as TRENDvent and a double line patient system for CARAT II

All systems have 22 mm cone connectors.

Order number:	
Tube system without integrated pressure measuring tube:	00007875
Tube system with integrated pressure measuring tube:	00007116
Single line patient circuit for ventilation:	00007967
Double line patient circuit for ventilation:	00007968



Tube system without integrated pressure measuring tube



Tube system with integrated pressure measuring tube



Single line patient circuit for ventilation




Double line patient circuit for ventilation

Accessories CPAP

Filter

Original filters assure a proper function and a long life of the CPAP device.


Order number:	
TREND series:	00007870
VECTOR et series:	00002029
TREND II series:	00002058
point:	00002110
CARAT:	00004880



TREND



VECTOR et



TREND II
TRENDvent



point



CARAT

Humidifiers

Humidifiers are intended to humidifying the therapy air flow which is supplied to the patient.

Order number:	
AquaTREND III-H5 for TREND-Series:	00007990
AquaDROP et for VECTOR et-Series:	00002015
AquaTREND uni for TREND II-Series:	00002073
Aquapoint for point:	00007849



AquaTREND III-H5



AquaDROP et



AquaTREND uni



Aquapoint

Transport- and Function Cases

The cases are specially developed for the specific device series. The cases protect the device optimal on transport. They offer ample space for the accessories like therapy tubes, masks or the power supply.

Order number:	
TREND series:	00007870
VECTOR et series:	00002029
TREND II series:	00002058
point:	00002110
CARAT:	00004880



TREND



point



VECTOR et



TREND II

For the ventilator series CARAT, two cases are offered.

Order number:	
Transport case CARAT I and II:	00004875
Function case CARAT I and II:	00004879
TRENDvent:	on request



CARAT I und II



TRENDvent

Accessories point

powerpack point

The powerpackpoint consists two batteries. Fully loaded the therapy device can operate for approx. 8 hours (without humidifier).

Order number: 00007846



filtersystem point

The ultrafilter is more efficient than the standard filter, because it is bigger and filters the respiratory therapy air twice.

Order number: 00007847

